



### Limited Warranty, Warranty Exclusions, Limit of Liability and RMA Request Instructions:

Before using, users shall determine the suitability of the product for their intended use, and users assume all risk and liability whatsoever in connection therewith. See the Menda Terms and Conditions - Mendapump.com/TermsAndConditions.aspx

# **MENDA**

#### Lasting-Touch Pump:

- One touch operation frees operator's hands for use
- Spill resistant design and swing-lid reduces fumes and waste
- Recommended for one-handed dispensing of liquids, solvents, and chemicals
- Durable, light-weight construction
- Compatible with any 45-400 threaded neck container
- Dispenses approximately 0.1cc of liquid with each pump

Pump Material	Polypropylene
Spring Material	Stainless-Steel
Locks for Travel	No
Neck Size / Thread	45/400
Lid-Type	Swing-lid

#### Pump Care:

Some types of fluids may leave a residue in the mechanism of the dispenser. Over time, this can reduce the performance of the pump. Occasional cleaning is recommended to maintain your unit to original efficiency.

Clean by flushing the mechanism thoroughly with warm soapy water. Flush the mechanism again with clean water to remove any soap residue. Do not boil or steam the dispenser. For long periods of inactivity, pumps should be flushed free of liquid and stored inverted in a dry place.

0	Specifications and procedures subject to change without notice. Unless otherwise noted, tolerance is ±10%.			
	LASTING-TOUCH, PUMP ONLY			
	3651 WALNUT AVE., CHINO, CA 91710 PHONE: (909) 627-2453 WEBSITE: <u>MendaPump.com</u>	DRAWING NUMBER 35600	DATE: April 2016	

## **Mouser Electronics**

Authorized Distributor

Click to View Pricing, Inventory, Delivery & Lifecycle Information:

Menda:

35600 35608 35601 35602 35603